



Standard operating procedure  
Passaging hiPS cells on mTeSR Plus + Vitronectin

Document code: iPS\_SOP\_0105.1

Page : 1/2

Author : M. Tok

Status :

Authorization : H.Mikkers

Authorization date: 12-10-2020

Print date : 24-4-2025

Revision date :

---

**Materials:**

- Vitronectin coated 6-well plate (iPS\_SOP\_0032)  
Date of coating:
- mTeSR Plus (iPS\_SOP\_0099)  
Prepared on:
- GCDR (SCT/07174)  
Lot: exp.:
- CellAdhere Dilution Buffer (SCT/07183)  
Lot: exp.:
- Cell scraper (Greiner/541-070)
- 10 µl filter tips (Corning/4807)
- 15 ml tube (Greiner/188 271)
- 1000 µl filter tips (Corning/4809)

***Important note:***

*Keep different lines in separate 6-well plates to avoid cross-contamination.*

*Be sure that all materials are at RT before starting the procedure.*

*Important:* *Check SCT mTESR Plus manual for details.*

**Method (based on one 6-well):**

- ☐ Remove differentiated parts of the colonies with a 10 µl pipet tip.
- ☐ Aspirate the medium and treat colonies with 1 ml GCDR for 2-5 min.  
*NOTE: Incubation time will differ per clone; colonies are ready when small wholes visible at 10x magnification start to appear within the colonies.*
- ☐ Remove GCDR and add 1 ml mTeSR Plus.
- ☐ Gently detach colonies by scraping with a cell scraper.
- ☐ Transfer the detached cell aggregates to a 15 mL tube with a P1000.
- ☐ Carefully breakup the colonies by pipetting the cell aggregate suspension up and down 2 - 3 times with a P1000. Avoid making single cells.



Standard operating procedure

Passaging hiPS cells on mTeSR Plus + Vitronectin

Document code: iPS\_SOP\_0105.1

Page : 2/2

Author : M. Tok

Status :

Authorization : H.Mikkers

Authorization date: 12-10-2020

Print date : 24-4-2025

Revision date :

- 
- ☐ Remove vitronectin from the vitronectin coated 6-well plate.
  - ☐ Add 2 ml mTeSR Plus medium to each well.
  - ☐ Add X  $\mu$ l of the cell aggregate suspension onto this new plate according the desired split ratio. Usual split ratio is between 1:20 – 1:100
  - ☐ Place the plate in a 37°C incubator. Move the plate in several quick, short, back-and-forth and side-to-side motions to disperse cells across the surface of the well.
  - ☐ Refresh every other day with 2 ml mTeSR Plus or 4 ml mTeSR Plus over the weekend.

**Note:** \_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_