

### **iPSCs differentiation into each of the three germ layers (endoderm, mesoderm and ectoderm)**

To determine if our cells are truly pluripotent, it is important to verify their ability to differentiate into each of the three germ layers (endoderm, mesoderm and ectoderm). In brief, approximately  $9 \times 10^3$  cells were plated per well of 8-well chamber slide, and grown overnight at 37 °C and 5% CO<sub>2</sub>, or until they reached 50% confluency. Then we used the Human Pluripotent Stem Cell Functional Identification Kit (SC027B, R&D Systems®) to differentiate into endoderm and ectoderm. In addition, to differentiate into mesoderm, we used the Differentiation Base Media Supplement (50X), that comes with this kit, with CHIR99021 supplement (STEMCELL™ Technologies), a WNT pathway activator that promotes direct mesoderm differentiation. We used the kit antibodies for immunocytochemical studies: Goat Anti-Human Otx2 Antigen Affinity-purified Polyclonal Antibody (ectoderm marker), Goat Anti-Human Brachyury Antigen Affinity-purified Polyclonal Antibody (mesoderm marker), and Goat Anti-Human SOX17 Antigen Affinity-purified Polyclonal Antibody (endoderm marker). Differentiated iPSCs were fixed in 4% paraformaldehyde (Sigma-Aldrich) for 10 minutes, incubated with phosphate-buffered saline solution with 0.1% Tween™ (PBST) and 1% BSA (Sigma-Aldrich) for 30 min and stained by standard immunofluorescence procedures. Cells were analyzed on DM400 M fluorescence microscope (Leica).