

ULEIi001-A

Differentiation potential – Ektoderm

Spontaneous differentiation of embryoid bodies.

Medium Essential 6 on Geltrex coated plates.

Antibodies:




- Neurofilament heavy Novus Bio NB300-135 rabbit polyclonal
- Secondary antibody Thermo Scientific goat anti-rabbit AF 488

Negative control – Gibco Episomal hiPSC line (A18945 Life Technologies)

Immunochemical stains were performed by MS Aniela Skrzypczyk




ULEIi001-A differentiation

Ektoderm – NEFH (neurofilament heavy)



ULEIi001-A differentiation

Ektoderm – NEFH (neurofilament heavy)



ULEIi001-A differentiation

Ektoderm – NEFH (neurofilament heavy)

