

Invoice number: SCTC2021-00079

# Certificate of Analysis 2021

Name principal investigator: Alex Garanto Cell line number: IPS18-00092 Project name: FFB

#### Table 1: Information on the reprogrammed cell line

Information cell line:	
Product description	Trilineage differentiation
Parental cell line	N/A*
Parental cell type	N/A*
Diagnosis	N/A*
Mutation	N/A*
Number of clones	1
Passage (P) of iPSCs reported at delivery	P11

\*N/A: Not Applicable

#### Table 2: Information on the characterization of the reprogrammed cell line

Test description:	Test method:	Test specification:	Result:
Activation of stem cell markers	qPCR	Upregulation of <i>SOX2, LIN28, NANOG, OCT4</i> in iPSCs compared with PBMCs	N/A*
Expression of stem cell markers	Immunocytochemistry	Expression of OCT4, NANOG, SSEA4, TRA-1-81	N/A*
Mycoplasma test	PCR	Negative	N/A*
Three lineage differentiation	Differentiation assay	Upregulaton of germlayer-specific markers	Pass

#### Three germ layer differentiation

IPS18-00092 clone 1 was differentiated into the endodermal, mesodermal and ectodermal germ layers. The RNA was isolated and the gene expression was checked by qPCR. Ct values are normalized with the housekeeping gene GUSB, set at 1. For each lineage two genes were assessed (Table 3). The differentiated cells were also stained for lineage-specific markers (Table 4).


#### Table 3: qPCR markers for three lineage differentiation

Lineage	Marker
Endoderm	FOXA2, SOX17
Mesoderm	Brachyury, HAND1
Ectoderm	PAX6, NCAM1

#### Table 4: ICC markers for three lineage differentiation

Lineage	Marker
Endoderm	SOX17
Mesoderm	NCAM1
Ectoderm	NESTIN


### Endoderm



# Upregulation of endodermal markers

Figure 1: Expression fold difference of endoderm-specific genes in differentiated cells, compared with undifferentiated iPSCs. *DNMT3B* was used as a reference for pluripotency.


## Mesoderm



Upregulation of mesodermal markers

Figure 2: Expression fold difference of mesoderm-specific genes in differentiated cells, compared with undifferentiated iPSCs. *DNMT3B* was used as a reference for pluripotency.

### Ectoderm



Upregulation of ectodermal markers

Figure 3: Expression fold difference of ectoderm-specific genes in differentiated cells, compared with undifferentiated iPSCs. *DNMT3B* was used as a reference for pluripotency.




Figure 4: Immunofluorescence staining of differentiated cells showing a positive signal of germlayer-specific markers.

Pass Fail Other:

Silvialbes

Silvia Albert, PhD Manager, Radboud Stem Cell Technology Center Date